

Fingertip Pulse Oximeter

MD300C2ES USER MANUAL

Ver1.0

General Description

Oxygen Saturation is a percentage of Oxyhemoglobin (HbO₂) capacity, compounded with oxygen, by all combinative hemoglobin (Hb) capacity in blood. In other words, it is consistency of Oxyhemoglobin in blood. It is a very important parameter for the Respiratory Circulation System. Many respiratory diseases can result in oxygen saturation being lowered in human blood. Additionally, the following factors can reduce oxygen saturation: Automatic regulation of organ dysfunction caused by Anesthesia, Intensive Postoperative Trauma, injuries caused by some medical examinations. That situation might result in light-headedness, asthenia, and vomiting. Therefore, it is very important to know the oxygen saturation of a patient so that doctors can find problems in a timely manner.


The fingertip pulse oximeter features low power consumption, convenient operation and portability. Place one fingertip into the photoelectric sensor for diagnosis and the pulse rate and oxygen saturation will appear on the display. It has been proven in clinical experiments that it also features high precision and repeatability.

Measurement Principle

Principle of the oximeter is as follows: A mathematical formula is established making use of Lambert Beer Law according to Spectrum Absorption Characteristics of Reductive hemoglobin (RHb) and Oxyhemoglobin (HbO₂) in red and near-infrared zones. Operation principle of the instrument: Photoelectric Oxyhemoglobin Inspection Technology is adopted in accordance with Capacity Pulse Scanning and Recording Technology, so that two beams of different wavelength of lights (660nm red and 905nm near infrared light) can be focused onto a human nail tip through a clamping finger-type sensor. A measured signal obtained by a photosensitive element, will be shown on the oximeter's display through process in electronic circuits and microprocessor.

Diagram of Operation Principle

1. Red and Infrared-ray Emission Tube
2. Red and Infrared-ray Recept Tube



Precautions For Use

1. Before use, carefully read the manual.
2. Operation of the fingertip pulse oximeter may be affected by the use of an electrosurgical unit (ESU).
3. The fingertip pulse oximeter must be able to measure the pulse properly to obtain an accurate SpO₂ measurement. Verify that nothing is hindering the pulse measurement before relying on the SpO₂ measurement.
4. Do not use the fingertip pulse oximeter in an MRI or CT environment.
5. Do not use the fingertip pulse oximeter in situations where alarms are required. The device has no alarms. It is not for continuous monitoring.
6. Do not use the fingertip pulse oximeter in an explosive atmosphere.
7. The fingertip pulse oximeter is intended only as an adjunct in patient assessment. It must be used in conjunction with other methods of assessing clinical signs and symptoms. In order to ensure correct sensor alignment and skin integrity, the maximum application time at a single site for our device should be less than half an hour.
8. Do not sterilize the device using autoclaving, ethylene oxide sterilizing, or immersing the device in liquid. The device is not intended for sterilization.
9. Follow local ordinances and recycling instructions regarding disposal or recycling of the device and device components, including batteries.
10. This equipment complies with IEC 60601-1-2:2007 for electromagnetic compatibility for medical electrical equipment and/or systems. However, because of the proliferation of radio-frequency transmitting equipment and other sources of electrical noise in healthcare and other environments, it is possible that high levels of such interference due to close proximity or strength of a source might disrupt the performance of this device.
11. Portable and mobile RF communications equipment can affect medical electrical equipment. This equipment is not intended for use during patient transport outside the healthcare facility.
12. This equipment should not be used adjacent to or stacked with other equipment.
13. It may be unsafe to:
 - use accessories, detachable parts and materials not described in the instructions for use
 - interconnect this equipment with other equipment not described in the instructions for use
 - disassemble, repair or modify the equipment.
14. These materials that contact with the patient's skin contain medical silicone and ABS plastic enclosure are all pass the ISO10993-5 Tests for invitro cytotoxicity and ISO10993-10 Tests for irritation and delayed-type hypersensitivity.

Contraindication

It is not for continuous monitoring.

Inaccurate measurements may be caused by

1. Significant levels of dysfunctional hemoglobin (such as carboxyl - hemoglobin or methemoglobin).
2. Intravascular dyes such as indocyanine green or methylene blue.
3. High ambient light. Shield the sensor area if necessary.
4. Excessive patient movement.
5. High-frequency electrosurgical interference and defibrillators.
6. Venous pulsations.
7. Placement of a sensor on an extremity with a blood pressure cuff, arterial catheter, or intravascular line.
8. The patient has hypotension, severe vasoconstriction, severe anemia, or hypothermia.
9. The patient is in cardiac arrest or is in shock.
10. Fingernail polish or false fingernails.
11. Weak pulse quality (low perfusion).
12. Low hemoglobin.

Product Properties


1. Simple to operate and convenient to carry.
2. Small volume, light weight and low power consumption.
3. High brightness OLED display SpO₂, PR, and Pulse bar; two directions display mode.
4. 2 display modes.
5. 2pcs AAA-size alkaline batteries; battery-low indicator.

Intended Use

The Fingertip Pulse Oximeter is only for sports and aviation use. It is ideal for use during sports activities, mountain climbing and piloting airplanes. It is not intended to diagnosis any medical condition or to be used in medical applications.


Operation Instructions

1. Install two AAA batteries according to the Battery Installation instructions.
 2. Place one of your fingers into the rubber opening of the pulse oximeter.
 3. Press the switch button one time on front panel to turn the pulse oximeter on.
 4. Keep your hands still for the reading. Do not shake your finger during the test. It is recommended that you do not move your body while taking a reading.
 5. Read the data from the display screen.
 6. After take out the finger, the measurement data displays in the screen for 5 seconds and then the pulse oximeter will power off automatically in 5 seconds.
 7. Press the button again to toggle between two display modes.
- After turning on the Oximeter, each time you press the power switch, the Oximeter will switch to another display mode. There are 2 display modes shown as follows:



Holding the power switch for longer than one second, will adjust the brightness of the oximeter. There are 10 levels of brightness. The default level is level four.

Front Panel




The pulse bar less than 30% indicates signal inadequacy and the displayed SpO₂ and pulse rate value is potentially incorrect.

Battery Installation

1. Install two AAA batteries into the battery compartment. Match the plus (+) and minus (-) signs in the compartment. If the polarities are not matched, damage may be caused to the oximeter.
2. Slide the battery door cover horizontally along the arrow shown as the picture.

Note:

- ◆ Please remove the batteries if the pulse oximeter will not be used for long periods of time.
- ◆ Please replace the battery when the low battery power indicator.



Using the Lanyard

1. Thread thinner end of the lanyard through the loop.

2. Thread thicker end of the lanyard through the threaded end before pulling it tightly.

Warnings!

- ◆ Keep the oximeter away from young children. Small items such as the battery door, battery, and lanyard are choking hazards.
 - ◆ Do not hang the lanyard from the device's electrical wire.
1. Please notice that the lanyard which is tied to the oximeter may cause strangulation due to excessive length.

Maintenance and Storage

1. Replace the batteries in a timely manner when low voltage lamp is lighted.
2. Clean surface of the fingertip oximeter before it is used in diagnosis for patients.
3. Remove the batteries if the oximeter is not operated for a long time.
4. It is best to store the product in -20°C ~ +55°C and ≤93% humidity.
5. Keep in a dry place. Extreme moisture may affect oximeter lifetime and may cause damage.
6. Dispose of battery properly; follow any applicable local battery disposal laws.

Cleaning the fingertip pulse oximeter

Please use medical alcohol to clean the silicone touching the finger inside of oximeter with a soft cloth dampened with 70% isopropyl alcohol. Also clean the being tested finger using alcohol before and after each test.

Do not pour or spray liquids onto the oximeter, and do not allow any liquid to enter any openings in the device. Allow the oximeter to dry thoroughly before reuse.

The use life of the device is five years when it is used for 15 measurements every day and 10 minutes per one measurement. Stop using and contact local service center if one of the following cases occurs:

- An error in the *Possible Problems and solutions* is displayed on screen.
- The oximeter cannot be powered on in any case and not the reasons of battery.
- There is a crack on the oximeter or damage on the display resulting readings cannot be identified; the spring is invalid; or the key is unresponsive or unavailable.

Specifications

1. Display Type

OLED display

2. SpO₂

Display range: 0%~100%

Measurement range: 70%~100%


Accuracy: 70%~100%: ±2%; 0%~69% no definition

Resolution: 1%

ARMS Value Analysis

Testing Range	Data Points	Bias	ARMS	Specification
70--<80	81	-0.70	1.82	±2
80--<90	71	-0.60	1.52	±2
90--100	84	-0.15	1.78	±2
70~100	236	-0.47	1.72	±2

Bland-Altman graphical plot of all sampled data points as below



Note: A functional tester cannot be used to assess the accuracy of a pulse oximeter monitor or sensor. Clinical testing is used to establish the SpO₂ accuracy. The measured arterial hemoglobin saturation value (SpO₂) of the sensors is compared to arterial hemoglobin oxygen (SaO₂) value, determined from blood samples with a laboratory CO-oximeter. The accuracy of the sensors in comparison to the CO-oximeter samples measured over the SpO₂ range of 70%~100%. Accuracy data is calculated using the root-mean-squared (Arms value) for all subjects, per ISO 9919:2005, Medical Electrical Equipment - Particular requirements for the basic safety and essential performance of pulse oximeter equipment for medical use

A functional tester is used to measure how accurately Fingertip Pulse Oximeter is reproducing the specified calibration curve and the PR accuracy.

The model of functional tester is Index2 FLUKE simulator and the version is 2.1.3.

3. Pulse Rate

Display range: 0bpm~254bpm

Measure range: 30bpm~235bpm

Accuracy: 30bpm~99bpm, ±2bpm; 100bpm~235bpm, ±2%

Resolution: 1bpm

4. LED Wavelengths

	Wavelength	Radiant Power
RED	660 ± 3nm	3.2mw
IR	905 ± 10nm	2.4mw

NOTE: The information about wavelength range can be especially useful to clinicians.

5. Power Requirements

Two AAA alkaline Batteries

Power consumption: Less than 30mA

Battery Life: Two AAA 1.5V, 600mAh alkaline batteries could be continuously operated as long as 30 hours.

6. Environment Requirements

Operation Temperature: 5°C ~ 40°C


Storage/ Transport Temperature: -20°C ~ +55°C

Ambient Humidity: ≤80% no condensation in operation; ≤93% no condensation in storage/transport

Atmosphere pressure: 86kPa ~ 106kPa

7. Equipment data update period

As shown in the following figure. Data update period of slower average is 8s.



8. Classification

According to the type of protection against electric shock: INTERNALLY POWERED EQUIPMENT;

According to the degree of protection against electric shock: TYPE BF APPLIED PART, (applied part: the rubber hole of the device);

According to the degree of protection against ingress of water: IPX1

According to the mode of operation: CONTINUOUS OPERATION

Declaration

Guidance and Manufacturer's declaration – electromagnetic emissions-For all EQUIPMENT and SYSTEMS			
Guidance and Manufacturer's declaration - electromagnetic emission			
The Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of Pulse Oximeter should assure that it is used in such an environment.			
EMission test	Compliance	Electromagnetic Environment – guidance	
RF emissions CISPR 11	Group 1	The Pulse Oximeter uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.	
RF emissions CISPR 11	Class B	The pulse Oximeter is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.	
Harmonic emissions	Not Applicable		
IEC 61000-3-2			
Voltage fluctuations/ flicker emissions	Not Applicable		
IEC 61000-3-3			


Guidance and Manufacturer's declaration – electromagnetic immunity-For all EQUIPMENT and SYSTEMS			
Guidance and Manufacturer's declaration - electromagnetic immunity			
The Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of the Pulse Oximeter should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance Level	Electromagnetic Environment – guidance
Electrostatic Discharge (ESD)	+/- 6kV contact +/- 8kV air	+/- 6kV contact +/- 8kV air	Floors should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%.
IEC 61000-4-2			

Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be at levels characteristics of a typical location in a typical commercial or hospital environment.
---	------	------	--

Guidance and Manufacturer's declaration – electromagnetic immunity-For all EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING

Guidance and Manufacturer's declaration - electromagnetic immunity

The Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of the Pulse Oximeter should assure that it is used in such an environment.

Immunity test	IEC test level	60601 test level	Compliance Level	Electromagnetic Environment – guidance
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	3 V/m	<p>Portable and mobile RF communications equipment should be used no closer to any part of the Pulse Oximeter, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance</p> $d=1.2\sqrt{P}$ <p>80 MHz to 800 MHz: $d=2.3\sqrt{P}$ 800 MHz to 2.5 GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection structures, objects and people.

^a Field strengths from fixed transmitters, such as base station for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Pulse Oximeter should be observed to verify normal operation. If abnormal performance is observed, additional measurements may be necessary, such as reorienting of the relocating the Pulse Oximeter.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m

Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEMS - For all EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING

Recommended separation distances between portable and mobile RF communications equipment and Pulse Oximeter

The *Pulse Oximeter* is intended for use in electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the *Pulse Oximeter* can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the *Pulse Oximeter* as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)	
	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d=1.2\sqrt{P}$	$d=2.3\sqrt{P}$
0.01	0.1167	0.2334
0.1	0.3689	0.7378
1	1.1667	2.3334
10	3.6893	7.3786
100	11.6667	23.3334

For transmitters rated at a maximum output power not listed above, the recommended separation distance in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.
















NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Possible Problems and Solutions

Problems	Possible reason	Solution
SpO ₂ or PR cannot be shown normally	1. Finger is not inserted correctly 2. Patient's SpO ₂ value is too low to be measured	1. Retry by inserting the finger 2. There is excessive illumination 3. Try some more times. If you can make sure no problem exist in the product, please go to a hospital timely for exact diagnosis.
SpO ₂ or PR is shown unstably	1. Finger might not be inserted deep enough. 2. Excessive patient movement	1. Retry by inserting the finger 2. Be calmness
The oximeter cannot be powered on	1. No battery or low power of battery 2. Batteries might be installed incorrectly 3. The oximeter might be damaged	1. Please replace batteries 2. Please reinstall the batteries 3. Please contact with local customer service centre
Indication lamps are suddenly off	1. The product is automatically powered off when no signal is detected longer than 8 seconds 2. The battery power is too low to work	1. Normal 2. Replace the batteries
"Error3" or "Error4" is displayed on screen	1. Err 3 means the red emission LED is damaged 2. Err 4 means the infra-red emission LED is damaged	1. Please contact with local customer service centre 2. Please contact with local customer service centre
Error 6	Err 6 means the screen is failure	Please contact with local customer service centre
"Error7" is displayed on screen	Err 7 means all the emission LED or reception diode is damaged.	Please contact with local customer service centre

Symbol Definitions

Symbol	Definition	Symbol	Definition
	Type BF applied part.		Low power indication
	Caution		No SpO ₂ Alarm
	Protected against dripping water.		Serial No.
	Oxygen saturation		Manufacturer's information
	Pulse rate (BPM)		Date of Manufacture
	Storage temperature and relative humidity		European union approval
	Authorized representative in the European community		Follow instruction for use
	Waste electrical and electronic equipment		

Box Contents

- Fingertip pulse oximeter
- One lanyard
- Two AAA batteries
- One instruction manual
- One carrying pouch

Applicable Models

MD300C2ES MD300C2ES-6

LIMITED WARRANTY

1. Limited Warranty

This Limited Warranty is valid only if you purchased the product from a ChoiceMMed America Corporation ("ChoiceMMed") authorized Reseller.

Customers can obtain the Limited Warranty to the applicable warranty period specified for different ChoiceMMed's products:

The warranty period of equipment is two years. The warranty period of accessories is six months.

1.1 ChoiceMMed warrants the software, mechanical and electronic components of this product to be free of defects in material and workmanship if used under normal operating conditions.

If the product displays any defects within the specified warranty period and that defect is not excluded under clause 4, ChoiceMMed will either repair or replace the product using new or reconditioned product or parts. In the event that ChoiceMMed decides to replace the entire product, this Limited Warranty shall apply to the replacement product for a new specified warranty period from the date of replacing of the original product.

1.2 When a product or part is exchanged, any replacement item becomes your property and the replaced item becomes ChoiceMMed's property. Parts provided by us in fulfillment of its warranty obligation must be used in products for which warranty service is claimed.

1.3 To obtain warranty service, you must deliver the product, freight collect, in either its original packaging or packaging providing an equal degree of protection, to the address specified by ChoiceMMed. In accordance with applicable law, ChoiceMMed may require that you furnish proof of purchase details and/or comply with registration requirements before receiving warranty service. It is your responsibility to backup any data, software, or other materials you may have stored or preserved on the product. It is likely that such data, software, or other materials will be lost or reformatted during service, and ChoiceMMed will not be responsible for any such damage or loss.

2. Online Registration

Please do remember to register your new ChoiceMMed equipment, which can quickly be completed by visiting www.choicemmedamerica.com/register. We would ask that you read the terms and conditions of our Limited Warranty carefully. If the product purchase is registered, this will help you obtain our full warranty service at a convenient and efficient way.

3. Warranty Service; Return Authorization Number

If you are not able to solve your issue, you need to contact us to request a Return Authorization (RA) Number. All inquiries must be accompanied by a description of the problem, the RA number, and a copy of the original sales receipt.

RETURN ADDRESS: ChoiceMMed America Corporation C/O Keystone Industrial Park, 2558 Pearl Buck Rd, Suite 8A, Bristol, PA 19007

3.1 No returns can be accepted without a valid RA number.

3.2 Purchaser should not have to pay for freight if there is a warranty issue or a defect issue with the product

4. Warranty Limitations and Exclusions

4.1 This Limited Warranty does not cover consumable parts. These include, but are not limited to: display face, batteries, carrying case, lanyard, etc.

4.2 This Limited Warranty does not cover the product if it has been electronically or mechanically modified in any way provided that any modifications were not done by ChoiceMMed or their authorized personnel. If the product needs to be modified or adapted in order to comply with applicable local technical or safety standards in any country which is not the country for which the product was originally developed and manufactured, this modification / adaptation shall not be considered a defect in materials or workmanship. This Limited Warranty does not cover any such modification / adaptation, regardless of whether it was carried out professionally or not. Under the terms of this Limited Warranty, ChoiceMMed shall not be held responsible for any cost resulting from such a modification / adaptation provided that any modifications were not done by ChoiceMMed or their authorized personnel

4.3 Product damage or defects caused by the following conditions are not covered by this Limited Warranty:

4.3.1 improper handling, neglect or failure to operate the unit in compliance with the instructions given in user or service manuals

4.3.2 connection or operation of the unit in any way that does not comply with the technical or safety regulations applicable in the country where the product is used.

4.3.3 any Act of God or Nature (such accident, fire, flood, etc.) or any other condition that is beyond the control of ChoiceMMed.

4.4 Any repair or opening of the unit carried out by unauthorized personnel (including the user) will void the Limited Warranty completely.

4.5 Products which do not meet the terms of this Limited Warranty will be repaired exclusively at the buyer's expense. ChoiceMMed or its authorized service center will inform the buyer of any such circumstance and request a written order to repair the product. In the event that the buyer fails to submit a written repair order within 6 weeks of notification, ChoiceMMed reserves the right to dispose of the product in an appropriate, environmentally friendly manner.

5. Warranty Transferability

This Limited Warranty is extended exclusively to the original buyer (customer of authorized Reseller) and is not transferable to anyone who may subsequently purchase this product.

No other person (distributor, dealer, fulfiller, retailer etc.) shall be entitled to give any warranty promise on behalf of ChoiceMMed.

6. Other Warranty Rights and National Law

6.1 This Limited Warranty does not exclude or limit the buyer's statutory rights as a consumer in any way.

6.2 The Limited Warranty regulations mentioned herein are applicable unless they constitute an infringement of applicable statutory local laws.

6.3 This Limited Warranty does not detract from any statutory seller's obligations in regard to any lack of conformity of the product and any hidden defect.

7. Amendment

Warranty service conditions are subject to change with the approval of both parties.

Contact Information

Phone: 215 874 0458

Email: service@choicemmedamerica.com

Distributed by:

ChoiceMMed America Corporation
Bristol, PA 19007

US Customer Service Site

ChoiceMMed America Corporation
C/O Keystone Industrial Park
2558 Pearl Buck Rd: Suite 8A
Bristol, PA 19007
Phone: (215)874-0458
www.choicemmedamerica.com

ChoiceMMed America Corporation

©2015 All RIGHTS RESERVED

Revised date: December 15, 2015